



SURGICAL MASKS

1

WASH HANDS BEFORE AND AFTER
REMOVAL OF MASK

2

PHYSICALLY INSPECT MASK
BEFORE USE

3

ENSURE PROPER FIT - STARTS AT
NOSE BRIDGE AND ENDS DOWN
PAST THE MOUTH

4

AVOID TOUCHING INSIDE OF MASK

5

FOLD ENTRANCES SHOULD BE
FACING DOWNWARD

6

PLACE LEFT ELASTIC BEHIND LEFT
EAR

7

PLACE RIGHT ELASTIC BEHIND
RIGHT EAR

8

STORE MASK IN DESIGNATED AREA
(BROWN PAPER BAG) OR UNDER
CONTROLLED CONDITIONS (UV
LIGHT)

9

ENCOURAGE PRACTICAL HYGIENE
REUSE MASK

10

DISCARD APPROPRIATELY WHEN
USE HAS EXPIRED (24 HOURS)

EMAIL: MHATA@WESTERNU.EDU OR
KEVIN.MANCILLA@WESTERNU.EDU