



N-95 MASKS

N-95 masks should be prioritized for high risk patient contact

- 1 WASH HANDS BEFORE AND AFTER REMOVAL OF MASK
- 2 PHYSICALLY INSPECT MASK BEFORE USE
- 3 ENSURE PROPER FIT - NEEDS TO COVER NOSE AND MOUTH
- 4 PLACE LOWER STRAP BEHIND NECK
- 5 PLACE TOP STRAP ON BACK OF HEAD
- 6 AVOID TOUCHING INSIDE OF MASK
- 7 IF MASK CONTAINS METAL NOSE BAND, BEND METAL TO SHAPE AROUND NOSE WITHOUT PINCHING AIRWAY
- 8 STORE MASK IN DESIGNATED AREA (BROWN PAPER BAG) OR UNDER CONTROLLED CONDITIONS (UV LIGHT)
- 9 ENCOURAGE PRACTICAL HYGIENE REUSE OF N-95 MASKS - SHOULD NOT EXCEED 5 USES
- 10 DISCARD APPROPRIATELY WHEN USE HAS EXPIRED

REFERENCES: PLEASE VISIT
[HTTPS://MANAGEDCAREPHARMACYCA.COM](https://managedcarepharmacyca.com)



HOW TO CONSERVE AND REUSE MASKS

- 1 ENCOURAGE REUSE BY AVOIDING CONTACT WITH INTERIOR OF MASK
- 2 IF MULTIPLE N-95 MASKS: ROTATE MASKS AND DECONTAMINATE IN BETWEEN ROTATIONS
- 3 PRIORITIZE USE FOR HIGH RISK PATIENTS AND STAFF
- 4 STATION **ONE** STAFF MEMBER AT FRONT OF STORE WITH APPROPRIATE PPE TO SCREEN SYMPTOMATIC PATIENTS AND EMPLOYEES
- 5 PRIORITIZE SITUATIONS WHEN STAFF IS HAVING PATIENT CONTACT (I.E. VACCINATIONS)
- 6 IF USING N-95: SOME DATA SHOWS PLACING A SURGICAL MASK ON TOP MAY CONSERVE LONGEVITY OF THE N-95
- 7 AFTER USE: STORE N-95 MASKS IN BROWN PAPER BAG OR IN CONTROLLED UV LIGHT SETTING THAT IS CONTAMINATION FREE
- 8 CONSIDER CLOTH MASK IF NONE OTHER ARE AVAILABLE OR IF YOU MUST CONSERVE N-95 MASKS

REFERENCES: PLEASE VISIT [HTTPS://MANAGEDCAREPHARMACYCA.COM](https://managedcarepharmacyca.com)
EMAIL: MHATA@WESTERNU.EDU OR KEVIN.MANCILLA@WESTERNU.EDU



SURGICAL MASKS

- 1 WASH HANDS BEFORE AND AFTER REMOVAL OF MASK
- 2 PHYSICALLY INSPECT MASK BEFORE USE
- 3 ENSURE PROPER FIT - STARTS AT NOSE BRIDGE AND ENDS DOWN PAST THE MOUTH
- 4 AVOID TOUCHING INSIDE OF MASK
- 5 FOLD ENTRANCES SHOULD BE FACING DOWNWARD
- 6 PLACE LEFT ELASTIC BEHIND LEFT EAR
- 7 PLACE RIGHT ELASTIC BEHIND RIGHT EAR
- 8 STORE MASK IN DESIGNATED AREA (BROWN PAPER BAG) OR UNDER CONTROLLED CONDITIONS (UV LIGHT)
- 9 ENCOURAGE PRACTICAL HYGIENE REUSE MASK
- 10 DISCARD APPROPRIATELY WHEN USE HAS EXPIRED (24 HOURS)

EMAIL: MHATA@WESTERNU.EDU OR
KEVIN.MANCILLA@WESTERNU.EDU