



N-95 MASKS

N-95 masks should be prioritized for high risk patient contact

1

WASH HANDS BEFORE AND AFTER REMOVAL OF MASK

2

PHYSICALLY INSPECT MASK BEFORE USE

3

ENSURE PROPER FIT - NEEDS TO COVER NOSE AND MOUTH

4

PLACE LOWER STRAP BEHIND NECK

5

PLACE TOP STRAP ON BACK OF HEAD

6

AVOID TOUCHING INSIDE OF MASK

7

IF MASK CONTAINS METAL NOSE BAND, BEND METAL TO SHAPE AROUND NOSE WITHOUT PINCHING AIRWAY

8

STORE MASK IN DESIGNATED AREA (BROWN PAPER BAG) OR UNDER CONTROLLED CONDITIONS (UV LIGHT)

9

ENCOURAGE PRACTICAL HYGIENE REUSE OF N-95 MASKS - SHOULD NOT EXCEED 5 USES

10

DISCARD APPROPRIATELY WHEN USE HAS EXPIRED

REFERENCES: PLEASE VISIT

[HTTPS://MANAGEDCAREPHARMACYCA.COM](https://managedcarepharmacyca.com)